



Mounting frame insert Flex M, petridishes

- for petridishes d=88 mm, 65 mm, 54 mm, 36 mm
- for insertion in Mounting frame Flex M





Mounting frame insert Flex M, slides and chambers

- for object slides 76x26 mm
- for chambers (ibidi® µ-Slide™, Nunc® Lab-Tek™/Lab-Tek™ II)
- for insertion in Mounting frame Flex M




Mounting frame Flex M


- for microwell plates, flasks and petridishes d=95 mm
- for the mounting of Mounting frame inserts Flex M



Specimen stage 232x230 with hardcoat anodized surface *



Object guide M 130x85



Aqua Stop II
(optional in addition)
431716-0000-000

Illuminator TL LED 10 W
(included with the stand)
423004-9040-000

Illumination System RGB-UV, equipped with 4 solid state LEDs

- Red (625 nm)
- Green (565 nm)
- Blue (470 nm)
- UV (385 nm)

(included with the stand)

Objective LD A-Plan 5x/0.15 Ph1 M27 *

Objective LD A-Plan 10x/0.25 Ph1 M27 *

Objective LD A-Plan 20x/0.35 Ph1 M27 *

Objective LD A-Plan 40x/0.55 Ph1 M27 *

(further objectives by choice)

Filter Set 90 HE LED with Reflector Module FL EC P&C *
(further Filter Sets by choice)

Microscope Axiocvert 5 digital; integrated color camera
431030-9210-000

Microscope Axiocvert 5 digital; integrated mono camera
431030-9190-000

Light shield to block ambient light

* Specimen stage, objectives and Filter Set 90 HE LED with Reflector Module FL EC P&C are already pre-installed on the microscope

These components are included in the scope of delivery:

- External power adaptor
- Country-specific power cable
- USB 3.0 Type C data and 60 W fast charging cable for Apple® iPad Pro™ or Microsoft® Surface™



Dust protection set
434308-9010-000
(included with the stand)


(Apple and Microsoft are registered trade marks)

The following components are not included in the scope of delivery and must be procured separately:

Microsoft® Surface Pro 7™ (i5-1035G4 or i7-1065G7 processor) 12.3",
or: Apple® iPad Pro™ 4th Gen 12.9",
or: Apple® iPad Pro™ 5th Gen 12.9"

Keyboard

Mouse



Wenzel Mikroskoptechnik GmbH

Seit 1996

Ihr kompetenter Fachhändler für

[ZEISS-Mikroskope](#)

Telefon [03641-828080](tel:03641-828080)

eMail verkauf@mikroskoptechnik.de

im Netz www.mikroskoptechnik.de